

## Zeeman-effect

The Zeeman-effect contains all effects that can be observed related to spectral lines whose emitters are exposed to a magnetic field. It is differentiated between the "normal" and the "anomalous" Zeeman-effect. The normal Zeeman-effect emerges at pure orbital angular momentum  $L$  ( $S=0$ ), every state splits into  $2L+1$  equidistant levels. At the anomalous Zeeman-effect the spin- and the orbit-momentum couple to the total momentum  $J=L+S$  ( $S \neq 0$ ) causing a splitting into  $2J+1$  levels.

At strong external magnetic fields the Paschen-Back effect occurs whereby the L-S-coupling becomes negligible compared to the energy of the single magnetic moments in the field.

In the experiment the shift from  $3^3D$  to  $2^3P$  ( $\lambda = 587.6$  nm) in Helium in a variable magnetic field is analyzed.

### I. Required knowledge

Multiplet systems, nomenclature of terms (Lit. 1, 2, 3)

Zeeman effect (normal and anomalous) and Paschen-Back effect (Lit. 1, 2, 3)

Selection rules and polarization (Lit. 1, 2, 3)

Energy level schemes for the visible He-spectrum (Lit. 1)

Arrangement of the Fabry-Perot spectrometer (dispersion range, resolving power, interference order, finesse) (Lit. 4)

Digital camera (CMOS and CCD array, creation of colors) (Lit. 5, 6)

### II. Literature

1. H. Haken, Atomic and Quantum Physics, Springer-Verlag 1987  
H.C. Wolf:
2. W. Demtröder: Atoms, molecules and photons, Springer 2006
3. R. Eisberg Quantum physics of atoms, molecules, solids, nuclei and particles, 1985  
R. Resnick:
4. W. Demtröder: Laser spectroscopy, Springer 2003
5. <http://micro.magnet.fsu.edu/primer/digitalimaging/concepts/ccdanatomy.html>
6. <http://www.microscopyu.com/articles/digitalimaging/ccdintro.html>

### III. Experimental task

The splitting of a spectral line of Helium ( $3^3D-2^3P$ ,  $\lambda = 587.6$  nm) in a magnetic field is investigated at different field strengths with the help of a Fabry-Perot interferometer. The ring-shaped interference pattern is recorded by a digital camera.

1. Calibrate the magnetic field with the Hall probe (remanence of the magnet?).
2. Adjust the optical assembly of the transversal observation so that you can see a sharp, high-contrast and good illuminated ring system.
3. Vary the exposure time, shutter and ISO-Values for optimal results.
4. Take images of the ring system at different magnetic field strength (increment 50 mA – 100 mA, up to 1.0 T).
5. Adjust the optical assembly of the longitudinal observation so that you can see sharp, high-contrast and good illuminated sections of the ring system.
6. Take images by rotating the polarizer (increment  $5^\circ$ ).

### IV. Experimental analysis

Transversal observation:

Plot the squared of the ring diameters against the magnetic field  $B$ .

Determine Bohr's magneton  $\mu_B$  under the assumption, that Paschen-Back effect prevails.

Longitudinal observation:

Create a Video-clip from the images and explain what can be seen.

### V. Necessary data

Wavelength:  $\lambda = 587.6$  nm

Etalon distance:  $d = 9.462$  mm

Interference filter:  $\lambda = 589$  nm, HW = 23 nm


Pixel dimensions:  $7.81 \mu\text{m} \times 7.81 \mu\text{m}$

### VI. Attention

1. **Caution:** High voltage up to 4.9 kV (4.0 mA) at the lamp! Do not change the adjustment of the power supply and do not touch the set-up while the lamp is switched on.
2. **Caution:** The current for the magnet coils has to be lower than 1.0 A for permanent use, 2.0 A for short times.
3. Do not turn bolts and nuts at the Fabry-Perot-Etalon!  
The polarization of the turnable polarizer is **not** indicated by the mark on the mounting but by the red line near the screws (important for the optimal adjustment).

## VII. Experimental set-up:

### Transversal:



### Longitudinal:

